

XN-1000™

Change the way you
perform hematology



Cell analysis in hematology that's one step ahead

The XN-1000 is designed to be a step ahead of hematology analyzers of the past and change the way you perform hematology in the modern clinical laboratory. This compact system offers proven and innovative fluorescent flow cytometry technology that delivers a highly comprehensive CBC and differential.


Sysmex Automated Hematology Analyzers Offer Expanded Clinical Value

NRBCs and 6-part differential standard on all models:

- NRBC testing on every CBC; reduces the need for reflex testing or performing automated/manual differentials.
- WBC count does not require correction in presence of NRBCs.
- Automated reportable immature granulocyte count (IG#, IG%) with every automated differential reduces manual review rates.

Fluorescent Platelet Technology¹:

- Platelet-specific stain greatly reduces interferences from RBC or WBC fragments.
- Can be set up as a reflex test² with onboard rules or Caresphere™ Workflow Solution (WS).



Fluorescent Reticulocyte Technology:

- Fluorescent flow cytometry generates reportable RET# and RET%, Reticulocyte Hemoglobin equivalent (RET-He) and Immature Retic Fraction (IRF).

Qualitative Flagging of Immature WBCs:

- The addition of the WPC¹ channel, along with enhanced algorithms, provides better detection of abnormal lymphocytes and blast cells.

Body Fluid Mode¹:

- Provides WBC (WBC-BF), RBC (RBC-BF), Total Nucleated Cell count (TC-BF), and a two-part WBC differential.

Low WBC Analysis Mode:

- Low WBC mode reduces vote-outs and manual reviews.
- Can be set up as a reflex test² with onboard rules or Caresphere WS.


Pre-Dilute Mode:

- Allows testing to be performed on samples of whole blood that are diluted with CELLPACK™ DCL.



*1 Channel not available on all model configurations

*2 Automatic reflex sample analysis is not available on all module configurations









Specifications

Principles and Technology	Fluorescent flow cytometry using a semi-conductor laser and hydrodynamic focusing in dedicated channels
Parameters	<p>Whole Blood: WBC; RBC; HGB; HCT; MCV; MCH; MCHC; PLT (PLT-I, PLT-F¹); NEUT#, NEUT%; IG#, IG%; LYMPH#, LYMPH%; MONO#, MONO%; EO#, EO%; BASO#, BASO%; NRBC#, NRBC%; RDW-SD; RDW-CV; MPV; RET#, RET%; IRF, RET-He; IPF#¹, IPF%¹</p> <p>Body Fluid¹: TC-BF, WBC-BF, RBC-BF PMN #, PMN%, MN #, MN%</p>
Throughput	Whole Blood Mode: up to 100 CBC+DIFF samples/hour per module Body Fluid Mode ¹ : up to 40 samples/hour per module
Sample Aspiration Volume	Whole Blood Mode: 88 µL Low WBC Mode: 88 µL Pre-dilute Mode: 70 µL Body Fluid Mode ¹ : 88 µL
Linearity	WBC: 0.00 – 440.00 X 10 ³ /µL RBC: 0.00 – 8.60 X 10 ⁶ /µL PLT: 0 – 5,000 X 10 ³ /µL WBC-BF ¹ : 0.000 – 10.000 X 10 ³ /µL RBC-BF ¹ : 0.000 – 5.000 X 10 ⁶ /µL
Quality Control	<ul style="list-style-type: none"> XN CHECK™: complete tri-level QC product for all CBC, Diff, PLT, and RET parameters XN CHECK™ BF: bi-level QC product for Body Fluid channel Comprehensive QC files including current and new lot feature 99 total files including 5 XbarM files Insight™ web-based Interlaboratory Quality Assessment Program (IQAP) BeyondCareSM Quality Monitor (BCQM) provides real-time peer group analytics and advanced troubleshooting assistance in a simplified platform

¹Channel not available on all model configurations

Specifications

	XN-1000 R	XN-1000 BR	XN-1000 PR	XN-1000 BPR	XN-1000-010
					
Model ID	XN-1000-100-R-01	XN-1000-100-BR-01	XN-1000-100-PR-01	XN-1000-100-BPR	XN-1000-010
Components					
XN-10 or XN-20	XN-10	XN-10	XN-10	XN-10	XN-20
Automated reflex/rerun	0	0	0	Included	Included
Reticulocyte License (R)	Included	Included	Included	Included	Included
Body Fluid License (B)	0	Included	0	Included	Included
Fluorescent Platelet (P)	0	0	Included	Included	Included
WPC Channel	0	0	0	0	Included
Floor-standing wagon option	Available				
Analyzer Dimensions	20.5" (W) X 26.6" (D) X 33.0" (H); 150 LBS			25.4" (W) X 29.7" (D) X 33.7" (H); 172 LBS	
					
Analyzer on Wagon Dimensions	25.4" (W) X 36.8" (D) X 61.5" (H); 350 LBS				

